
PediCap[®]

Pediatric End-Tidal CO₂ Detector

English	2
Français	6
Deutsch	10
Nederlands	14
Español	18
Italiano	22
Português	26
Svenska	30

PediCap[®]

Pediatric End-Tidal CO₂ Detector

 **Rx ONLY**

Directions for Use

Description

The Pedi-Cap[®] detector, when connected between an endotracheal tube and a breathing device, detects approximate ranges of end-tidal CO₂ by color comparison in patients weighing 1–15 kg (2.2–33 lb). The detector may be used during patient transport or in any location where intubations are performed. It may be used for up to two hours.

Indications

Use to assist verification of tube placement during endotracheal or nasotracheal intubation.

Use on intubated patients to detect approximate ranges of end-tidal CO₂ when clinically significant.

Contraindications

- Not to be used for detection of hypercarbia.
- Not to be used to detect main stem bronchial intubation.
- Not to be used during mouth-to-tube ventilation.
- Should not be used to detect oropharyngeal tube placement. Standard clinical assessment must be used.

WARNING:

- 1. Before using the Pedi-Cap detector, read complete Directions for Use.**
- 2. Do not use this device on patients with body weight greater than 15 kg due to potential increased airway resistance.**
- 3. The Pedi-Cap detector is not suitable for interpretation by individuals with blue-yellow color blindness.**

Instructions for Use

1. Match the initial color of the indicator to the purple color labeled CHECK around the detector window as shown in **A**.
If the purple color of the indicator is not the same color or darker than the area marked CHECK, DO NOT USE.
2. Insert the endotracheal tube. Inflate the cuff, if equipped.
3. Firmly attach the Pedi-Cap detector to the endotracheal tube; then attach the breathing device, as shown in **B**.
4. When a Heat and Moisture Exchanger (HME) is used, note:

- Connect the HME to the endotracheal tube
 - Connect the Pedi-Cap detector to the HME
 - Connect the breathing device to the Pedi-Cap detector
5. Ventilate the patient with SIX BREATHS of moderate tidal volume (may be done slowly or quickly). Compare the indicator color in the window on full-end expiration to ranges printed on the detector cover. Interpret results according to the decision tree that follows.

If Pedi-Cap detector results are not conclusive, the endotracheal tube should be immediately reinserted unless correct anatomic location can be confirmed with certainty by other means.

WARNING: INTERPRETING RESULTS BEFORE CONFIRMING SIX BREATH CYCLES CAN YIELD FALSE RESULTS.

Gastric distention with air prior to attempted intubation may introduce CO₂ levels as high as 4.5% into the Pedi-Cap detector if the endotracheal tube is misplaced in the esophagus. Initial Pedi-Cap detector color (yellow) may be interpreted as a false positive if read before delivery of six breaths.



WARNING: THE PEDI-CAP DETECTOR WILL RESPOND TO CO₂ IN THE PRESENCE OF A MAIN STEM BRONCHIAL INTUBATION. STANDARD CLINICAL ASSESSMENT SHOULD BE USED TO CONFIRM PROPER POSITION OF THE ENDOTRACHEAL TUBE WITHIN THE TRACHEA.

CAUTION

 **STORE AT OR BELOW 24°C (75°F). STORAGE FOR EXTENDED PERIODS AT TEMPERATURES ABOVE 24°C (75°F) MAY REDUCE SHELF LIFE.**

CAUTIONS:

1. Inspect the Pedi-Cap detector before attaching it to the endotracheal tube. Do not use if any physical damage is observed.
2. Large airway leaks may decrease tidal volume delivered to the detector which can limit the accuracy of the Pedi-Cap detector.
3. Do not use device in the presence of trichlorethylene or chloroform anesthetic; chemical interactions will affect device accuracy.

4. **Do not subject device to extreme humidity, such as a heated humidifier or nebulizer; excessive humidity will affect device accuracy.**
5. **Reflux of gastric contents, mucus, edema fluid, or intratracheal epinephrine into the Pedi-Cap detector can yield persistent patchy yellow or white discoloration which does not vary with the respiratory cycle. Contamination of this type may also increase airway resistance and affect ventilation. Discard device if this occurs.**
6. **In cardiac arrest, re-establishment of cardiac output and pulmonary perfusion by adequate CPR is necessary to increase end-tidal CO₂ to levels detectable by the Pedi-Cap detector (above 0.5%).**
7. **The Pedi-Cap detector is not a substitute for observation of the patient. This device must not be relied upon as the sole indicator of resuscitation performance.**

Mechanical Specifications

Internal Volume:	3 cc
Resistance to Flow:	2.5 cm H ₂ O ± 0.5 cm at 10 L/min flow
Weight:	Less than 5 g
Connector Ports:	Patient end: 18 mm O.D. / 15 mm I.D. Circuit end: 15 mm O.D. / 5 mm I.D.

Nellcor Puritan Bennett and Pedi-Cap are trademarks of Nellcor Puritan Bennett Incorporated.

Covered by one or more of the following U.S. Patents and foreign equivalents: 4,728,499; 4,879,999; 5,166,075; and 5,179,002.

© 2006 Nellcor Puritan Bennett Incorporated. All rights reserved.

Interpretation

Color Range A

End-Tidal CO₂ level approximately 0.03% to < 0.5%. < 4 mmHg, < 0.5 kPA


Color Range B

End-Tidal CO₂ level approximately 0.5% to < 2%. 4 to < 15 mmHg, 0.5–2 kPA


Color Range C

End-Tidal CO₂ level approximately 2% to 5%. 15 to 38 mmHg, 2–5 kPA

Pedi-Cap Interpretation in the Presence of Adequate Pulmonary Blood Flow/Spontaneous Heartbeat:

Pedi-Cap Interpretation in the Presence of Poor Pulmonary Blood Flow/Cardiac Arrest:

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© 2006 Nellcor Puritan Bennett Inc. Alla rättigheter förbehålles.
Tyco Healthcare UK LTD